

Antibody Specification Sheet

Mouse IgG1 - Purified Isotype Control

Catalog reference: AC1-P
Size: 0.1 mg



Antibody Information:

Antigen: unknown.

Ig Class: mIgG1 .

Form: purified.

Specificity: unknown.

Antibody Source: mono-clonal antibody from a mouse myeloma.

Production: *in vitro* cell culture.

Purification: Protein A agarose chromatography.

Purity: ≥95%

Formulation: provided as a 0.2 µm sterile-filtered solution in Ca⁺⁺ & Mg⁺⁺ free Dulbecco's

PBS with 0.05% sodium azide.

Concentration: 1 mg/ml.

Recommendations: For use as an isotype control.

Storage conditions: store undiluted at 4 °C or in aliquots at ≤ -20 °C.

Applications:

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Antibody Solutions

1135 San Antonio Road, Suite D., Palo Alto, CA 94303, USA
Tel: 650-938-4300, FAX: 650-938-4390
E-mail: solns@antibody.com