

## Antibody Specification Sheet

### Mouse IgG2b - Purified Isotype Control

---

Catalog reference: AC3-P

Size: 100 µg



#### Antibody Information:

**Antigen:** unknown.

**Ig Class:** mIgG2b .

**Form:** purified.

**Specificity:** unknown.

**Antibody Source:** mono-clonal antibody from a mouse myeloma.

**Production:** *in vitro* cell culture.

**Purification:** Protein A agarose chromatography.

**Purity:** ≥95%

**Formulation:** provided as a 0.2 µm sterile-filtered solution in Ca<sup>++</sup> & Mg<sup>++</sup> free Dulbecco's

PBS with 0.05% sodium azide.

**Concentration:** 1 mg/ml.

**Recommendations:** For use as an isotype control.

**Storage conditions:** store undiluted at 4 °C or in aliquots at ≤ -20 °C.

#### Applications:

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

### Antibody Solutions

1135 San Antonio Road, Suite D., Palo Alto, CA 94303, USA

Tel: 650-938-4300, FAX: 650-938-4390

E-mail: [solns@antibody.com](mailto:solns@antibody.com)